

Silipos polymer gel technology provides a unique dermal-delivery system that gradually diffuses USP medical grade mineral oil to impart deep moisturisation for soft and supple skin. Silipos gels have been developed for numerous Medical and Beauty Care applications, and can be customized with a wide range of additives to meet specific consumer needs.

As the leading, vertically integrated manufacturer of patented polymer gel technologies in the world, many leading global corporations rely on Silipos for their product development needs. All gel development, including the incorporation of additives, is conducted in-house, as well as the base gel manufacturing and conversion to finished products.

Another unique feature of Silipos gels are the "visco-elastic" properties which allow the gel to move synergistically with the skin helping to prevent abrasion and discomfort. And of course, all Silipos gels are non-toxic, latex free and hypoallergenic.

Silipos' orthopedic and foot care products are made with medical grade mineral oil gels that easily conform to the body ensuring a comfortable fit that cushions and protects. When applied directly to the skin, Silipos gels have the added benefit of skin moisturising and conditioning. The orthopedic and foot care range offers easy to use solutions for a variety of ailments ranging from common foot discomforts to post traumatic injuries. Whatever your needs, Silipos has a product that can help improve the quality of life for your patients.

## Sizing Chart - Tubing, Pads & Caps

All sizing is approximate

CODE	DESCRIPTION	FLAT WIDTH	LENGTH	RELAXED	STRETCHED
10055	Digital Pad Mini	18mm	22mm	51mm	70mm
10065	Digital Pad S/M	19mm	25mm	57mm	89mm
10075	Digital Pad L/XL	29mm	32mm	76mm	114mm
1105	Digital Pads on a Strip Mini	18mm	610mm	51mm	70mm
1106	Digital Pads on a Strip S/M	19mm	610mm	57mm	89mm
1107	Digital Pads on a Strip L/XL	29mm	610mm	76mm	114mm
1100	DigiStrip Narrow	16mm	610mm	45mm	70mm
1101	DigiStrip Medium	19mm	610mm	51mm	89mm
1102	DigiStrip Wide	25mm	610mm	64mm	101mm
10855	Gel Tubing Mesh Narrow	19mm	152mm	45mm	70mm
10845	Gel Tubing Mesh Wide	76mm	254mm	152mm	254mm
10865	Gel Tubing Ribbed Narrow	19mm	152mm	45mm	67mm
10875	Gel Tubing Ribbed Wide	29mm	152mm	64mm	108mm
10195	Digital Cap Ribbed Mini	16mm	64mm	38mm	64mm
10205	Digital Cap Ribbed S/M	19mm	64mm	48mm	76mm
10206	Digital Cap Ribbed L/XL	25mm	64mm	64mm	95mm
10235	Digital Cap Mesh	19mm	64mm	41mm	64mm
6511	Gel Toe Separators Small	38mmL x 23mmW x 4mmH			
6512	Gel Toe Separators Large	54mmL x 23mmW x 4mmH			
6508	Gel Toe Spreaders Small	24mmL x 12mmW x 20mmH			
6509	Gel Toe Spreaders Medium	28mmL x 16mmW x 25mmH			
6510	Gel Toe Spreaders Large	31mmL x 16mmW x 32mmH			