

# SlimLine™ Cast Boot

## Low Profile Traditional Cast Boot



The DARCO SlimLine™ Cast Boot can be adapted to fiberglass casts, bulky compression bandages and a variety of woundcare applications, making it one of the most versatile cast shoes on the market.

**DARCO**



## SlimLine™ Features and Benefits:

- › **Square Toe Design** acts as a protective bumper and gives the patient more room to accommodate larger casts and bandages.
- › **Forefoot Closure** can be adjusted to accommodate casts or bandages and can be switched between left/right closure.
- › **Higher Ankle Strap** adds stability.
- › **Rocker Sole** allows for a smooth roll throughout the gait pattern.

## SlimLine™ Indications:

- › Use with casts and compression bandaging of the foot and ankle
- › Surgery of the midfoot/forefoot

Item Number	Description	Size	Color	Women's Size	Men's Size
SLQ0B	SlimLine Cast Boot	X-Small	Black	3.5 - 5	
SLQ1B	SlimLine Cast Boot	Small	Black	5.5 - 8	
SLQ2B	SlimLine Cast Boot	Medium	Black	8.5 - 10	6 - 8
SLQ3B	SlimLine Cast Boot	Large	Black	10.5 - 13	8.5 - 11
SLQ4B	SlimLine Cast Boot	X-Large	Black	13.5+	11.5 - 14

Available through:



Headquarters:

DARCO International | 810 Memorial Blvd. | Huntington, WV 25701 | USA

Toll Free voice 800-999-8866 | Fax 877-201-7082 | Voice 304-522-4883 | Fax 304-522-0037

info@darcointernational.com | www.darcointernational.com

*As with all DARCO products, a portion of the profits from the sale of these items will be donated to encourage and support continuing education in our profession.*