





Smart-V-Link Software

Manufactured in:
Japan



- Allows for creation of detailed patient reports in pdf format - ideal for referrals
- Fully auto-gain controlled waveform by computer
- ABI, TBI and arterial blood flow velocity
- Real-time vascular waveform display
- Data storage for future reference
- Standardised testing modules for easy operation and documentation
- Includes specially designed USB interface cable to link Doppler to PC
- Requires Windows 7/8.1/10 (laptop not included)

DOP/SMLINK

Log-in to www.briggatemedical.com for pricing

Minimum System Requirements:

Compatible with the following Hadeco models:

- ES100V3
- Smartdop 45
- Bidop 7
- Smartdop 30EX

Operating System: Windows 10/8.1/7

CPU: Celeron 1 GHz or more

Display: 800 x 600 dots+, 256 colours

Memory: 1Gb+ RAM, 40MB of HDD space

USB: USB 2.0 required as Software is supplied on USB Flash Drive
Doppler connection via Unique USB Cable (supplied)

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com