



Smartdop 30EX

BI-DIRECTIONAL DOPPLER

Warranty:
1 year

Manufactured in:
Japan




- LCD Real-time Waveform Display
- Includes: 8MHz Probe, Leg Cuff & In-built Printer
- Automatic Cuff Inflation and Deflation
- 30 Waveform Memory
- Dual Probe Connection
- Consolidates ABI (ankle-brachial index) and TBI (toe-brachial index) systolic pressure calculations
- Automatic pulse volume read-outs for segmental leg study
- Automatic Venous Reflex readouts

DOP-SMA30

Log-in to www.briggatemedical.com for pricing

Dimensions:		Compatible with:	
Width	25.2cm	- PPG Probe	
Depth	19.2cm	- Smartlink Software	
Height	9.1cm	- Digital Cuff	
Weight	1.6kg	- Temperature Probe	

ADD SMART-V-LINK SOFTWARE WHEN PURCHASING THE HADECO SMARTDOP 30EX & SAVE \$75.00+GST



briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggatemedical.com

www.briggatemedical.com