



Smartdop 45

BI-DIRECTIONAL DOPPLER

Warranty:
1 year

Manufactured in:
Japan




- Bidirectional Hand-Held Doppler
- Backlit LCD Screen displays realtime waveforms and numerical data
- 8MHz Probe
- Convenient Probe Button freezes data for notation
- In-built Printer for Waveform Documentation
- 30 Waveform Memory
- Ni-MH Battery with Charger
- Speaker Output: 250 mW +
- Includes Carry Bag and Headphones

DOP-SMA45

Log-in to www.briggatemedical.com for pricing

Dimensions:		Compatible with:	
Width	9.2cm	-	PPG Probe
Depth	21cm	-	Smartlink Software
Height	6cm	-	Temperature Probe
Weight	440g		

ADD SMART-V-LINK SOFTWARE WHEN PURCHASING THE HADECO SMARTDOP 45 & SAVE \$75.00+GST



briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro) | 1800 33 4142 (Free Call) | briggate@briggate.com.au | www.briggatemedical.com

© Copyright Briggate Medical Company 2020 - No parts of this document either written or images may be reproduced without permission.