

Suda Mobility Dom. Case



Multiple external pockets to neatly house consumables

Tri-Wheel design makes ascending and descending stairs and gutters a breeze



Specifications:

- **Work (Drill) Unit**
 - External Dimensions: 35cmW x 25cmH x 35cmD
 - Internal Dimensions: 33cmW x 19cmH x 25cmD
 - Weight: 4kg (empty)
- **Supply Unit**
 - 35cmW x 48cmH x 35cmD
 - Weight: 5.8kg (empty)
- **Trolley**
 - 52cmW x 107cmH x 35cmD
 - Weight: 2.7kg (empty)

2 x Plastic Drawer units slide snugly into partitions within the lower "Supply Unit" case

Trolley frame is collapsable to help facilitate transport between patients

VIDEO DEMO:



Suda Upright Podiatry Case with Trolley

Includes:

- Upper "Work Unit" Case is removable - designed to hold Podiatry Drill
- Lower "Supply Unit" Case is removable - contains 6 x Drawers to hold loose items yet still has room for additional consumables and instruments
- Adjustable L.E.D. Work Light
- Cases sit upon the "Mobility Trolley" with Tri-Wheel configuration to make ascending and descending stairs a breeze

CASE-SUDA1 Suda Upright Trolley Case

Log-in to www.briggatemedical.com for pricing



Removable/adjustable L.E.D Work Light

"Work Unit" Case has a slot in the rear and a pocket to house the power lead to help reduce set-up times.

Multiple internal pockets to neatly house gloves, dressings and other consumables

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com