

SÜDA vac S

briggatedrills
SALES • SERVICE • SUPPORT

PERFORMANCE
GUARANTEE

contact us for further information



Suda Vac 'S' Podiatry Drill with Dust Extraction

- 25,000 rpm Handpiece
- LCD Display
- Handpiece On/Off Switch
- Clockwise/Counter-clockwise Rotation
- Brushless Vacuum Motor Technology
- Touch Button Controls
- 3 Programmable Memory Buttons for Suction Level, Burr Speed & Burr Direction
- Customisable Handpiece Holder

- Height 12cm
- Width 28.5cm (includes h/p holder)
- Depth 23cm
- Weight 3.8kg
- Handpiece Weight 110g

SUDA-VACS

Suda Vac 'S' Podiatry Drill
with Dust Extraction

Log-in to www.briggatemedical.com for pricing

Suction Performance:

210 L/min

"Silent Suction"

Warranty:

2 years

Manufactured in:

Germany



briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com