

Swann-Morton®



Founded in Sheffield during 1932, Swann-Morton have become a world leader in the manufacture of surgical blades, scalpels and handles. It is a name respected globally for quality, precision, consistency and reliability and recognised for combining traditional blade making skills with the latest technology.

Microbiological

Microbiological tests are performed on packaging materials and various other materials used in the manufacturing process. Bioburden estimates are randomly carried out on the materials and finished product from the Packing Department. No post sterilisation tests are performed as the product release is parametric. Testing for pyrogenicity is not required.

Latex

All Swann-Morton products are Latex Free.

Surgical Blades

No.3 and No.4 fitting surgical blades are manufactured in accordance with BS 2982 and BS EN 27740. Boxes are marked with lot number and expiry dates for traceability. All Sterile Product boxes are bar coded. Sterile blades are sterilised by Gamma radiation to a minimum dose of 25 kGy.

Surgical Handles No.3 & No.4 fitting surgical handles are manufactured in accordance with BS 2982 and BS EN 27740. All other handles are manufactured to in-house standards. Boxes are marked with lot number and expiry dates for traceability.

Disposable Scalpels

Packaging is marked with lot number and expiry dates for traceability. Boxes are bar coded. Sterile disposable scalpels are sterilised by Gamma radiation to a minimum dose of 25 kGy

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com