

Domiciliary Drills



Saeshin Thumb Domiciliary Podiatry Drill

- Variable Speed (to 30,000 rpm) Micro-motor Handpiece
- Forward/Reverse Burr Control
- High torque
- Quiet operation
- Minimal Handpiece vibration
- Rechargeable Battery
(up to 10 hours operation time)

■ Height	3.6cm
■ Width	6.5cm
■ Depth	12.9cm
■ Weight	388g
■ Handpiece Weight	172g

THUMB

Saeshin Thumb
Domiciliary Podiatry Drill

Log-in to www.briggatemedical.com for pricing



Warranty:
12 Months

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com