

## Hear what podiatrists say about the TOURNI-COT™

“The Tourni-Cot saves me time and hassle. I don't bother with other improvised tourniquet methods anymore. It's just that simple.”

-- Dr. Stavros Alexopoulos  
Chicago, IL

“The Tourni-Cot fully exsanguinates fluid from the surgical field so the phenol is more effective. As a result of hemostasis, the Tourni-Cot has decreased the recurrence of nail growth following nail removal procedures. ”

-- Dr. Michael David  
Grand Rapids, MI

“I have been using the Tourni-Cot for over 10 years. It is sterile, easy and effective. It provides even and safe compression. Utilizing a product like the Tourni-Cot as opposed to a penrose or a rubber band communicates to our patients and staff that patient safety comes first in our practice”

-- Dr. Neal Mozen  
Southfield, MI



MAR-MED CO.

The **Tourni-Cot™** makes it safe and easy to create a bloodless field for procedures on digits. Hand surgeons, **podiatrists**, emergency rooms and operating theaters around the world have used the **Tourni-Cot™** for over 20 years.

**briggate**medicalcompany

Supporting Podiatrists with integrity and trust since 1987

**Overview:**

The Tourni-Cot™ is the exsanguinating digit tourniquet. Tourni-Cots make it both safe and easy to create a bloodless field during digit procedures. Tourni-Cot™ rings are simple to apply and exsanguinate by exerting pressure as they roll proximally onto a digit. Once in place the ring applies reliable pressure to allow for a visible field. The Tourni-Cot™ has been safely used by hand surgeons, podiatrists, in emergency rooms and operating theaters around the world for over 20 years.

**Patient Safety:**

Use of the Tourni-Cot™ is a reliable practice for achieving the intended healthcare outcome. The United Kingdom's National Health Service agreed after performing a study on the risk of other improvised methods. The Tourni-Cot's four color-coded sizes make it easy to select the correct ring for your procedure. And with bright coloration and an attached warning tag, the Tourni-Cot™ is almost impossible to forget on a digit. These factors advance patient safety and reduce medical errors, while being cost effective and easy to adopt into an improved workflow.

**Doctor Preference:**

Doctors prefer using the Tourni-Cot™ because it is easy to use and saves time. Using the Tourni-Cot communicates that patient safety comes first, unlike the use of a penrose drain, surgical glove, or other improvised device. A hand surgeon specially designed the Tourni-Cot™ due to the burden and risk of other methods. For over 20 years doctors have agreed that the Tourni-Cot™ improves care by reliably creating a bloodless field, and helping procedures go smoothly.

**Product Specifics:**

The Tourni-Cot™ is available in small, medium, large, and extra-large sizes. All Tourni-Cots come sterile and individually sealed in a Tyvek® pouch. A warning label referencing instructions for use is attached to each ring. Tourni-Cots are made from an elasticized silicone polymer, making them completely latex free, non-allergenic, and extremely durable. Upon completion the Tourni-Cot™ is removed by severing the ring at its notch.

**Quality Assurance:**

Mar-Med Co. maintains a quality system to ensure the provision of safe and effective medical devices. The Tourni-Cot™ is manufactured in the USA under strict compliance with the ISO 13485:2003 and FDA standards, and is CE marked for compliance with European regulations.



**Safety First**

"CE marked digital tourniquets which are labeled and/or brightly coloured should be used, in accordance with manufacturers' instructions. Surgical gloves should not be used as tourniquets."

"While the cost of CE marked tourniquets is higher than a pair of sterile surgical gloves, the cost of litigation claims [resulting from the use of surgical gloves] can run to five or six figures..."


From: *Reducing risks of tourniquets left on after finger and toe surgery - a study by the UK's National Health Service*  
- December 2009



**briggate**medicalcompany

Supporting Podiatrists with integrity and trust since 1987

**The Exsanguinating Digit Tourniquet**




- Easy to use
- Saves time
- Effective
- Sterile
- Low risk
- Latex-free



The **Tourni-Cot™** makes it **safe and easy** to create a **bloodless field** for procedures on digits. Hand **surgeons, podiatrists**, emergency rooms and operating theaters **around the world** have used the Tourni-Cot™ for over **20 years**.

METHOD	EASE OF USE	INSTANT BLOODLESS FIELD?	RISK OF PRESSURE INJURY	RISK OF FORGETTING?
Penrose Drain with Clamp	Low	No	High	Low
Rolled Surgical Glove	Low	No	Low	Yes
Rubber Band	Low	No	High	High
Pneumatic Tourniquet	Low	No	Low	Low
<b>Tourni-Cot by Mar-Med Co</b>	<b>High</b>	<b>Yes</b>	<b>Low</b>	<b>Low</b>



**Mar-Med Tourni-Cot Digit Tourniquet - Sterile**

- Simple, exsanguinating digit tourniquet
- Provides reliable and gentle pressure adequate to occlude the vessels
- LATEX-FREE
- Images shown are approx. actual size.

<b>TOURNS</b>	Small	Yellow
<b>TOURNM</b>	Medium	Orange
<b>TOURNL</b>	Large	Green
<b>TOURNXL</b>	Extra Large	Blue

Log-in to [www.briggatemedical.com](http://www.briggatemedical.com) for pricing.

**briggatemedical**company

Supporting Podiatrists with integrity and trust since 1987

# Tourni-Cot Large



07-2015

REF

TCL

LOT

07-12-12-A

QTY 1

Single use only. Do not reuse or resterilize.

STERILE EO



CE 0050

**MAR-MED CO.**

333 Fuller Ave. NE, Grand Rapids, MI 49503, USA  
www.marmedco.com

EC REP

Canonbury Products Ltd, 2 St. James Road,  
Brackley NN13 7XY, UK

LB-TCL Rev 00

SINGLE-USE ONLY.  
DO NOT REUSE OR  
RESTERILIZE.



THIS MEANS NO  
RE-PROCESSING

