

allcare

UltraPlus

ULTRASOUND



UltraPlus Therapeutic Ultrasound

Utilises ultrasound waves to induce a deep heat to localised tissue which increases relaxation and recovery to the site. The UltraPlus is a high quality, extremely functional, and visually eye-catching unit that provides the professional physician with an effective ultrasonic therapy. It's small size and light weight make it simple to transport and easy to use beyond its dedicated place of operation.

- 1 & 3 MHz detachable applicators
- Compact and easy to use
- 10 pre-set treatment programs that can be customised to suit user
- LCD display

ULTRAPLUS

UltraPlus Ultrasound

Log-in to www.briggatemedical.com for pricing

Specifications:

- Channel: 2 independent
- Intensity (max): 3W/cm² ± 20%
- Frequency: 1MHz/3MHz
- Range: 0.2w/cm - 3.0w/cm
- Treatment time: 1 to 10 minutes
- Power supply: Adaptor, input 110V-240V; 50/60Hz, output 15V/3A

Dimensions:

- Height 7.5cm
- Width 25cm
- Depth 18.5cm
- Weight 2.4kg

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com