

# Mobile Domiciliary Case

## Upright Domiciliary Case on Wheels

- Suitable for all Podiatry Drills available from Briggate Medical
- External Dimensions: 400mmW x 500mmH x 300mmD (980mmH with Handle fully extended)
- Drill Cavity Dimensions: 350mmW x 120mmH x 250mmD
- Weight: 7.5kg

### CASE-MOBIL

Large Upright Podiatry Domiciliary Case on Wheels

Log-in to [www.briggatemedical.com](http://www.briggatemedical.com) for pricing



Top compartment is suitable for storage of consumables or instruments

Drill cavity that fits all currently stocked Drills from Briggate Medical

4 Drawers that can be fastened shut to ensure contents are safe during transit

Retractable handle and wheels for ease of transport

Fastening clips and zip-up sides for safe and secure transit



**briggate**medicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

[briggate@briggate.com.au](mailto:briggate@briggate.com.au)

[www.briggatemedical.com](http://www.briggatemedical.com)