

walker's PROFESSIONAL GEL

Introducing Walker's GEL - an expanded selection of gel protection products designed to relieve common foot, ankle, hand and limb ailments. A variety of advanced gel formulations are utilised to provide the best possible outcomes for your patients - all of which are:



Ribbed Knit

A fabric weave that is typically stretchy & elastic, allowing it to expand & contract without distorting the shape. The ribbed effect makes the fabric feel a little thicker providing added cushioning.



Mesh Knit

A fabric weave defined by its lightweight feel & permeable texture. Unlike most types of fabric, which feature closely woven textures, mesh is woven loosely, resulting in thousands of tiny holes for high breath-ability.

SIZING CHART

CODE	DESCRIPTION	SIZE	DIMENSIONS
WG1006	Fabric Toe Spreader - Ribbed	Medium	25mmW x 35mmL
WG1007	Fabric Toe Spreader - Ribbed	Large	33mmW x 38mmL
WG1010	Digital Pads - Ribbed	Mini	20mmW x 22mmL
WG1011	Digital Pads - Ribbed	Small/Medium	23mmW x 25mmL
WG1012	Digital Pads - Ribbed	Large/X-Large	28mmW x 32mmL
WG1031	All Gel Toe Spacer w/Toe Loop	Medium	Ø 12mm
WG1032	All Gel Toe Spacer w/Toe Loop	Large	Ø 18mm
WG1037	Hammer Toe Prop	Small	16mmW x 29mmL x 20mmH Ø 12mm
WG1038	Hammer Toe Prop	Medium	19mmW x 32mmL x 23mmH Ø 13mm
WG1039	Hammer Toe Prop	Large	21mmW x 37mmL x 30mmH Ø 17mm
WG1042	Looped Toe Spacer	Small/Medium	11mmW x 35mmL x 15mmH Ø 10mm
WG1044	Looped Toe Spacer	Large/X-Large	35mmW x 45mmL x 22mmH Ø 17mm
WG1050	Digital Caps - Ribbed	Mini	22mmW x 64mmL

continued on next page...

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com

SIZING CHART

CODE	DESCRIPTION	SIZE	DIMENSIONS
WG1051	Digital Caps - Ribbed	Small/Medium	26mmW x 64mmL
WG1052	Digital Caps - Ribbed	Large/X-Large	36mmW x 67mmL
WG1056	Digital Caps - Mesh	Small/Medium	15mmW x 60mmL
WG1075	All-Gel Digital Caps	Small/Medium	22mmW x 48mmL
WG1076	All-Gel Digital Caps	Large/X-Large	30mmW x 51mmL
WG1100	All-Gel Toe Separators	Small	22mmW x 37mmL (2mm thickness)
WG1102	All-Gel Toe Separators	Large	23mmW x 53mmL (3mm thickness)
WG1125	All-Gel Toe Spreaders	Small	12mmW x 25mmL x 21mmH
WG1126	All-Gel Toe Spreaders	Medium	14mmW x 28mmL x 28mmH
WG1127	All-Gel Toe Spreaders	Large	17mmW x 28mmL x 30mmH
WG1152	Fabric Gel Strip - Ribbed		15mmW x 60cmL
WG1050	Fabric Gel Strip - Ribbed		20mmW x 60cmL
WG1051	Fabric Gel Strip - Ribbed		25mmW x 60cmL
WG1177	Gel Pads on a Strip - Ribbed	Mini	Ø 15mm Pad
WG1175	Gel Pads on a Strip - Ribbed	Small/Medium	Ø 20mm Pad
WG1176	Gel Pads on a Strip - Ribbed	Large/X-Large	Ø 25mm Pad
WG1200	Digital Tubing - Ribbed	Narrow	20mmW x 16cmL
WG1201	Digital Tubing - Ribbed	Wide	25mmW x 16cmL
WG1202	Digital Tubing - Mesh	Narrow	15mmW x 15cmL
WG1203	Fabric Gel Sleeve - Mesh	Fully-Coated	72mmW x 25cmL (flat, unstretched)
WG1209	Fabric Gel Sleeve - Mesh	Half-Coated	72mmW x 25cmL (flat, unstretched)
WG1276	Thin Forefoot Pad	One Size	107mmW x 111mmL
WG1290	Plantar Fasciitis Sleeves	One Size	72mmW x 110mmL (flat, unstretched)
WG1303/1306	Bunion Sleeve	Small/Medium	75mmW x 85mmL (flat, unstretched)
WG1304/1307	Bunion Sleeve	Large/X-Large	90mmW x 85mmL (flat, unstretched)
WG13112	Bunion Relief Sleeve	Small/Medium	75mmW x 85mmL (flat, unstretched)
WG13114	Bunion Relief Sleeve	Large/X-Large	90mmW x 85mmL (flat, unstretched)
WG1326S	Gel Metatarsal Pad	Small	85mmW x 40mmL (flat, unstretched)
WG1327S	Gel Metatarsal Pad	Medium	95mmW x 45mmL (flat, unstretched)
WG1342	Gel Forefoot Cover	Small/Medium	80mmW x 110mmL (flat, unstretched)

continued on next page...

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

SIZING CHART

CODE	DESCRIPTION	SIZE	DIMENSIONS
WG1343	Gel Forefoot Cover	Large/X-Large	100mmW x 130mmL (flat, unstretched)
WG1400	Gel Achilles Heel Pad	Small/Medium	76mmW x 190mmL (flat, unstretched)
WG1401	Gel Achilles Heel Pad	Large/X-Large	90mmW x 240mmL (flat, unstretched)
WG1405	Ankle Protection Sleeve	One Size	100mmW x 215cmL (flat, unstretched)
WG1412	Boot Protection Sleeve	Small/Medium	90mmW x 175cmL (flat, unstretched)
WG1413	Boot Protection Sleeve	Large/X-Large	100mmW x 220cmL (flat, unstretched)
WG1350	Gel Body Discs		Ø 102mm
WG1354	Gel Body Discs		Ø 63mm
WG1352	Gel Squares	Non-Adhesive	100mmW x 100mmL
WG1353	Gel Squares	Adhesive	100mmW x 100mmL
WG1360	Gel Dots	Self-Adhesive	Ø 25mm
WG2255/2275	Uni-Strap	Small/Medium	90mmW x 90cmL (flat, unstretched)
WG2265/2285	Uni-Strap	Large/X-Large	105mmW x 95cmL (flat, unstretched)
WG3004/3005	Gel Met Pad	Small/Medium	70mmW x 50cmL (Ø 15mm toe loop)
WG3006/3007	Gel Met Pad	Large/X-Large	85mmW x 60cmL (Ø 20mm toe loop)
WG3025/3026	All-Gel Toe Prop	Small	60mmL (Ø 17mm toe loop)
WG3027/3028	All-Gel Toe Prop	Medium	70mmL (Ø 17mm toe loop)
WG3029/3030	All-Gel Toe Prop	Large	75mmL (Ø 17mm toe loop)
WG5025S	Heel Cups	Small	40mmW x 65mmL
WG5026S	Heel Cups	Medium	45mmW x 85mmL
WG5027S	Heel Cups	Large	55mmW x 100mmL
WG5031S	Heel Cups	X-Large	60mmW x 115mmL
WG7050S	Full Length Insole	Small	235mmL
WG7051S	Full Length Insole	Medium	250mmL
WG7052S	Full Length Insole	Large	260mmL
WG8020/8021	Carpal Gel Sleeve	Small	75mmW x 165mmL (flat, unstretched)
WG8025/8026	Carpal Gel Sleeve	Large	90mmW x 215mmL (flat, unstretched)

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987