



Rapid-don

Flexi-fit

Specifications Sheet



Nitrile Examination Gloves

- Comfort, dexterity & tactility of Latex
- Safer & stronger than Latex
- Latex-free
- Powder-free
- Textured finger-tip
- Beaded Cuff
- Ambidextrous
- Colour: Light Blue

Product Dimensions:

Size	Palm Width (mm)	Length (mm)
X-Small	75 ± 5mm	240min
Small	85 ± 5mm	240min
Medium	95 ± 5mm	240min
Large	106 ± 5mm	240min
Extra Large	115 ± 5mm	240min

Product Thickness: Single wall thickness (mm)

	Specification
Cuff	0.05-0.08
Palm	0.06-0.09
Finger	0.08-0.14
Beading	0.5min

Physical Properties:

	Un-aged	Aged
Force at break (N)	6min	6min

Residual Powder: 2.0 mg/glove maximum

Quality Specification:

Inspection Criteria	Sampling Plan/AQL
Dimension	13 pcs (Median)
Force at Break	13 pcs (Median)
Water Tight Test	G1 - AQL 1.5
Visual Inspection - Major	G1 - AQL 2.5
Visual Inspection - Minor	G1 - AQL 4.0
Residual Powder	n = 5 (average)
Packaging Appearance	n = 5 (AND)
Glove Physical Counting	n = 5 (AND)



Incorporating Rapid-don technology

a unique polymer coating on the inner surface of the glove that assists in donning the glove, especially when hands are damp

Now supplied as a box of

300 GLOVES



CONFORMS TO EN 455

CONFORMS TO ASTM D6319

ISO 14001


ARTG# 197779

Re-Order Codes:

GLFLEX3RXS	X-Small
GLFLEX3RS	Small
GLFLEX3RM	Medium
GLFLEX3RL	Large
GLFLEX3RXL	Extra Large

Packaging:

300 Gloves per Box
10 Boxes of 300 per CARTON



Country of Manufacture: Malaysia

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com