



Flexi-fit nitrile

Specifications Sheet



Nitrile Examination Gloves - Powder Free

- Comfort, dexterity & tactility of Latex
- Safer & stronger than Latex
- Latex-free
- Powder-free
- Textured finger-tip
- Beaded Cuff
- Ambidextrous
- Colour: Purple



Product Dimensions:

Size	Palm Width (mm)	Length (mm)
X-Small	75 ± 5mm	240min
Small	85 ± 5mm	240min
Medium	95 ± 5mm	240min
Large	106 ± 5mm	240min
Extra Large	115 ± 5mm	240min

Re-Order Codes:

GLFLEXIXS	X-Small
GLFLEXIS	Small
GLFLEXIM	Medium
GLFLEXIL	Large
GLFLEXIXL	Extra Large

Packaging:

200 Gloves per Box
10 Boxes of 200 per CARTON

Country of Manufacture: Malaysia

Product Thickness: Single wall thickness (mm)

	Specification
Cuff	0.05-0.08
Palm	0.06-0.09
Finger	0.08-0.14
Beading	0.5min

Quality Specification:

Chemotherapy Drugs and Concentration (Tested for Resistance to Permeation by Chemotherapy Drugs as per ASTM D6978-05-Test Report PN 106336)

Chemotherapy Drugs and Concentration	Minimum Breakthrough Detection Time (minutes) (Specimen 1/2/3)
Carboplatin, 10 mg/ml (10,000 ppm)	>240 minute
Carmustine (BCNU), 3.3 mg/ml (3,300 ppm)	15.0 (15.0,15.1,15.2)
Cisplatin, 1.0 mg/ml (1,000 ppm)	>240 minute
Cyclophosphamide (Cytoxan), 20.0 mg/ml (20,000 ppm)	>240 minute
Dacarbazine (DTIC), 10.0 mg/ml (10,000 ppm)	>240 minute
Doxorubicin Hydrochloride, 2.0 mg/ml (2,000 ppm)	>240 minute
Etoposide (Toposar), 20.0 mg/ml (20,000 ppm)	>240 minute
Fluorouracil, 50.0 mg/ml (50,000 ppm)	>240 minute
Ifosfamide, 50.0 mg/ml (50,000 ppm)	>240 minute
Methotrexate, 25.0 mg/ml (25,000 ppm)	>240 minute
Mitomycin C, 0.5 mg/ml (500 ppm)	>240 minute
Mitoxantrone 2 mg/ml (2,000 ppm)	>240 minute
Paclitaxel (Taxol), 6.0 mg/ml (6,000 ppm)	>240 minute
Thiotepa, 10.0mg/ml (10,000 ppm)	30.1 (30.5,30.3,30.1)
Vincristine Sulfate, 1.0mg/ml (1,000 ppm)	>240 minute

Physical Properties:

	Un-aged	Aged
Force at break (N)	6min	6min

Residual Powder: 2.0 mg/glove maximum

Quality Specification:

Inspection Criteria	Sampling Plan/AQL
Dimension	13 pcs (Median)
Force at Break	13 pcs (Median)
Water Tight Test	G1 - AQL 1.5
Visual Inspection - Major	G1 - AQL 2.5
Visual Inspection - Minor	G1 - AQL 4.0
Residual Powder	n = 5 (average)
Packaging Appearance	n = 5 (AND)
Glove Physical Counting	n = 5 (AND)



Briggate Medical Company

A Division of Drewson International Pty Ltd
ABN: 92 006 698 559
23-25 Lakewood Blvd Braeside Victoria Australia 3195
briggate@briggate.com.au
www.briggatemedical.com

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com