



Rapid-don

Flexi-fit

Specifications Sheet



Nitrile Examination Gloves - Powder Free

- Comfort, dexterity & tactility of Latex
- Safer & stronger than Latex
- Latex-free
- Powder-free
- Textured finger-tip
- Beaded Cuff
- Ambidextrous
- Colour: Light Blue

Product Dimensions:

Size	Palm Width (mm)	Length (mm)
X-Small	75 ± 5mm	240min
Small	85 ± 5mm	240min
Medium	95 ± 5mm	240min
Large	106 ± 5mm	240min
Extra Large	115 ± 5mm	240min

Product Thickness: Single wall thickness (mm)

	Specification
Cuff	0.05-0.08
Palm	0.06-0.09
Finger	0.08-0.14
Beading	0.5min

Physical Properties:

	Un-aged	Aged
Force at break (N)	6min	6min

Residual Powder: 2.0 mg/glove maximum

Quality Specification:

Inspection Criteria	Sampling Plan/AQL
Dimension	13 pcs (Median)
Force at Break	13 pcs (Median)
Water Tight Test	G1 - AQL 1.5
Visual Inspection - Major	G1 - AQL 2.5
Visual Inspection - Minor	G1 - AQL 4.0
Residual Powder	n = 5 (average)
Packaging Appearance	n = 5 (AND)
Glove Physical Counting	n = 5 (AND)

Incorporating Rapid-don technology

a unique polymer coating on the inner surface of the glove that assists in donning the glove, especially when hands are damp



CONFORMS TO EN 455

CONFORMS TO ASTM D6319

ISO 14001

ARTG# 197779

Re-Order Codes:

GLFLEXRXS	X-Small
GLFLEXRS	Small
GLFLEXRM	Medium
GLFLEXRL	Large
GLFLEXRXL	Extra Large

Packaging:

200 Gloves per Box
10 Boxes of 200 per CARTON

Country of Manufacture: Malaysia



Briggate Medical Company

A Division of Drewson International Pty Ltd
ABN: 92 006 698 559
23-25 Lakewood Blvd Braeside Victoria Australia 3195
briggate@briggate.com.au
www.briggatemedical.com

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com