



Lite-fit

Specifications Sheet



Nitrile Examination Gloves - Powder Free

- Comfort, dexterity & tactility of Latex
- Safer & stronger than Latex
- Latex-free
- Powder-free
- Textured finger-tip
- Beaded Cuff
- Ambidextrous
- Colour: White

Product Dimensions:

Size	Palm Width (mm)	Length (mm)
X-Small	75 ± 5mm	240min
Small	85 ± 5mm	240min
Medium	95 ± 5mm	240min
Large	106 ± 5mm	240min
Extra Large	115 ± 5mm	240min

Product Thickness: Single wall thickness (mm)

	Specification
Palm	≥0.05
Finger	≥0.08
Weight	3.2g

Physical Properties:

	Un-aged	Aged
Tensile Strength(MPa)	≥14	≥14
Ultimate Elongation (%)	≥550	≥400

Residual Powder: 2.0 mg/glove maximum

Quality Specification:

Inspection Criteria	Sampling Plan/AQL
Dimension	S2 - AQL 4.0
Physical Properties	S2 - AQL 4.0
Water Tight Test	G1 - AQL 1.0
Visual Inspection - Major	G1 - AQL 2.5
Visual Inspection - Minor	G1 - AQL 4.0
Residual Powder	n = 5 (average)
Packaging Labelling	n = 5 (AND)
Glove Physical Counting	n = 5 (AND)

LOW Dermatitis Potential

manufactured without the use of accelerator chemicals and sulphur resulting in lower dermatitis potential.



CONFORMS TO EN 455

CONFORMS TO ASTM D6319

ISO 14001

ARTG# 197779

Re-Order Codes:

GLLITEXS	X-Small
GLLITES	Small
GLLITEM	Medium
GLLITEL	Large
GLLITEXL	Extra Large

Packaging:

200 Gloves per Box
10 Boxes of 200 per CARTON

Country of Manufacture: Malaysia



Briggate Medical Company

A Division of Drewson International Pty Ltd
ABN: 92 006 698 559
23-25 Lakewood Blvd Braeside Victoria Australia 3195
briggate@briggate.com.au
www.briggatemedical.com

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com