

Water Distiller



Water Distiller

- A cost effective source of high quality distilled water for all your sterilisation needs
- Stainless steel tank
- Produces 1 litre per hour
- Includes Water Storage Cannister, Charcoal Filter and Chamber Cleaning Crystals
- 1 year limited warranty

Specifications:

- External Dimensions: 29cm x 29cm x 390cmH
- Chamber Dimensions: 18cm x 20cm
- Weight: 3.5kg
- Power: 750W

DISTILL2

Water Distiller 4 Litre

Log-in to www.briggatemedical.com for pricing

Water Distiller consumables are also available. Refer to our website for further info.

IMPORTANT:
PLEASE READ INSTRUCTION MANUAL BEFORE OPERATING THE
WATER DISTILLER FOR THE FIRST TIME.
FAILURE TO DO SO MAY RESULT IN MACHINE MALFUNCTION,
REDUCE THE LIFE-SPAN OF THE MACHINE AND/OR VOID WARRANTY.
TO PROLONG THE WORKING LIFE OF THE UNIT, WE SUGGEST THAT
AN EXTERNAL TIMING DEVICE (NOT INCLUDED) BE USED TO
AVOID THIS DEVICE BOILING DRY DURING OPERATION.

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com