

New at DARCO

# Therapeutic Socks in the TASO Function series | Toe Alignment Sock

TASO: uncomplicated & cost effective solution for easing big toe pain as a result of Hallux Valgus. The TASO Toe Alignment Socks are the latest addition to our series in corrective therapies here at DARCO.

TASO Function **Comfort**  
Has an integrated corrective tape



TASO Function **Control**  
Has an adjustable strap



The TASO Function Comfort and the TASO Function Control are therapeutic socks designed to counter the effects of Hallux Valgus. Whether it's a conservative approach to alleviating big toe discomfort or a post-operative alignment solution, the TASO series is the perfect answer. With two variations of therapy, you can choose which one is right for you.




Therapeutic Socks in the TASO Function series is a clever choice for any individual who wishes to proactively support the alignment of the toes. Whether you wear high-heels or have had Hallux Valgus surgery, the TASO Toe Alignment Socks are the right choice for you.

## TASO Function Comfort Hallux Valgus Therapeutic Socks



**OPTION NOT  
AVAILABLE  
IN AUSTRALIA**



When the big toe is flexible, an 8° shift on the MTP1 can be achieved on average

## TASO Function Control Hallux Valgus Therapeutic Socks with Adjustable Strap



Custom re-positioning through an integrated adjustable strap

## Conservative and Preventative Care

Hallux Valgus is purported to be the result of improper, repetitive, biomechanical malfunctions exerted on the foot.

The TASO Comfort has a stream-line, laminated strip built into the fabric of the sock along the big toe, gently correcting the position of the joint into a biomechanically accurate alignment. This contra balance can ease the effects of mistreatment (abuse) on the foot during the day.

With the laminated corrective strip on the outside of the sock, there's no risk of skin irritation or discomfort. It's easy to put on and the wearability is high with the TASO function socks. You can give your feet the therapy they need, without even feeling it – 24/7.

The TASO Comfort is sold by the pair and both socks have the laminated corrective strip for easy and uncomplicated conservative and preventative therapy. Both socks are of the TABI design; meaning there is a separate pocket for the big toe.

## Postoperative positioning, protecting Hallux Valgus corrective surgery yet versatile enough for conservative therapy as well

Postsurgical toe alignment is vital to fortify the results of surgery. The TASO Control can act as a supportive, transitional reinforcement tool by facilitating the training of the toes to maintain proper positioning. The TASO Control is a suitable choice following the bandage/splinting phase after surgery.

The principle is simple as it is effective: A strong thin strap is neatly sewn to the sock at the boarder of the big toe. The soft yet rigid strap is lead around the back of the foot and is secured to a small segment integrated into the opposite side of the sock. This allows the individual to regulate just how much correction is applied to the big toe and when.

Although the TASO Control is sold in pairs, the corrective strap is built into one sock (Left or Right) and the other sock is neutral; completely without the corrective function. Both socks are of the TABI design; meaning there is a separate pocket for the big toe.

The TASO Function series is the first course of therapy for any individual to pro-actively off-set the negative effects of compromising footwear or as a light weight toe positioning support after surgery. The TASO socks are uncomplicated, comfortable and easy to wear making them the perfect solution for corrective therapy on a daily basis.

The TASO Function series offers you preventative therapy as well as a protective solution fortifying surgical procedures.

## Advantages of TASO Toe Alignment Socks:

- > The TASO socks do not push the toe into position (creating unwanted pressure points) but influences the alignment through gently pulling the toe.
- > The induced alignment proactively counters the effects of the malposition and reduces discomfort while still allowing for free movement of the toes.
- > TASO Therapeutic Socks bring the big toe into proper alignment and maintain that position.
- > There is no altering or hindering of the biomechanical function and the joint can move freely.
- > With regular use, the socks can counteract the development and the progress of potential deformities.

## Features of TASO Toe Alignment Socks:

- > Both Toe Alignment Socks can be worn day or night.
- > Both Toe Alignment Socks fit comfortably in most shoes where there is ample room for the toes and if desired, can also be worn under dress socks.
- > The socks can be washed like any other dress sock (up to 60° inside out).
- > The corrective laminated strips are on the outside of the socks for comfortable yet continuous therapy.
- > Made of a soft comb-cotton fabric.
- > Made in Germany
- > Available in 6 sizes



## Item List TASO



SIZE	ITEM NO.		
TASO Function	Comfort	Control right	Control left
35.0 – 36.0	TFCT1	TFCL1R	TFCL1L
37.0 – 38.0	TFCT2	TFCL2R	TFCL2L
39.0 – 40.0	TFCT3	TFCL3R	TFCL3L
41.0 – 42.0	TFCT4	TFCL4R	TFCL4L
43.0 – 44.0	TFCT5	TFCL5R	TFCL5L
45.0 – 46.0	TFCT6	TFCL6R	TFCL6L

**Colour: Black**

**TASO Toe Alignment Socks are sold in pairs:**

- > Comfort: Both left & right have the laminated corrective strip.
- > Control: The adjustable strap is provided on either the left sock or the right sock (not both and the second sock in the pair is neutral).